

*Effortless Equipment Visibility,
Everywhere You Need It.*

Advanced Real-Time Location Systems—Designed for Healthcare, Built to Scale

Equipment delays. Misplaced assets. Wasted staff time. Unnecessary rentals; the cost to care is real. In a fast-paced clinical environment, knowing the real-time location of your medical equipment isn't a luxury, it's essential. With CenTrak's advanced Bluetooth® Low Energy (BLE) and cloud-based asset tracking solution, healthcare organizations gain the visibility needed to streamline workflows, reduce costs, and ensure the right equipment is always available when and where it's needed.

Real-time alerts. Room-level precision. Scalable peace of mind.

CenTrak's BLE-powered real-time location system delivers more than just a dot on a map, it delivers operational clarity. Designed for cost-effective deployment across entire healthcare campuses, this robust solution ensures reliable asset visibility anytime, anywhere.

Key Features That Maximize Efficiency

- Real-time location and status for high-value mobile medical equipment
- Integrations with CMMS and EHR systems
- Zonal-to-room-level accuracy to support workflow automation
- Power-efficient BLE tags with multi-year battery life
- Compact designs to accommodate a wide variety of asset types
- Customizable reporting to inform purchasing and rental decisions

Experience BLE That Delivers More

CenTrak's platform adapts to your facility's needs and asset management strategy. Here's what sets CenTrak apart:



Room Level Location Certainty - Patented location and calibration algorithms deliver precise location data to support critical workflows with confidence.



Low Latency - Location refreshes every 15–20 seconds (compared to 1+ minute with some BLE systems).



Floor-Level Accuracy - Avoid "floor hopping" with CenTrak's vertical discrimination technology.



Fast, Scalable Deployment - Install and configure in ~5 minutes per room; 800,000 sq. Ft. facility live in under one week.



Reliable Performance - Healthcare location experts with technology tailored to complex clinical settings for seamless, consistent communication and proven system reliability.



Native Access Point Integrations - Maximize existing investments with native support for Cisco, Aruba, and other BLE-enabled infrastructure, no extra hardware required.



Two-Way BLE Communication - Over-the-air updates for quick and efficient maintenance and upgrades.

At a Glance: CenTrak BLE Asset Tags



BLE-Only Mini Asset Tag



BLE-Only Rugged Asset Tag

BLE-Only Mini & Rugged Asset Tag
Technical Specifications

Feature	Details
Battery	1-5 years (depending upon form factor and usage), field-replaceable
Form Factors	Mini, Rugged (IP67-rated/Waterproof)
Mounting Options	Adhesive, zip tie, screw

Accelerate visibility and ROI with CenTrak's BLE Asset Tracking.

Connect with a CenTrak expert today to schedule a live demo and see how quickly your facility can go live.

BLE Infrastructure & Tags



BLE Gateway



BLE Repeater



BLE-Only Mini Asset Tag



BLE-Only Rugged Asset Tag



BLE MM Staff Badge

Contact Us: www.centrak.com | marketing@centrak.com | 800-515-2928