

CenTrak integrates with a broad ecosystem of healthcare technologies to extend the value of RTLS across clinical, operational, communication, safety, facilities, and asset management workflows. Through its flexible integration platform, CenTrak enables organizations to connect RTLS data with new or existing enterprise systems using industry-standard protocols and interfaces, helping accelerate implementations while preserving investments in existing technologies.

Common protocols and standards supported include:

- REST APIs
- HL7v2 (mllp)
- FHIR (json)
- HL7v3 (XML)

The integration partners listed represent examples of known CenTrak integrations and are not intended to be an exhaustive list. Additional integrations may be available through standard protocols, middleware platforms, APIs, or customer-specific implementation requirements.

EMR/EHR

Cerner	Meditech
Epic (Epic Toolbox Designation)	

Security / Access Control / CCTV

AMAG	Johnson Controls
Genetec	C-CURE 9000
	Lynx Systems

CMMS / Asset Management

Accruent / Connectiv	Medimizer
AIMS	Nuvolo
EQ2	PeopleSoft
Four Rivers (TMS)	

Patient Engagement / Smart Room

CareSpace	OneView
Epic MyChart	pCare
eVideon (Vibe)	Sonifi
GetWell Network	

Nurse Call / Clinical Communications

Amplion Clinical Communications	InTouch Health
Ascom	Jeron
Austco	Rauland AllTouch
Cisco Patient Connect	Rauland Enterprise
Connexall	Spok (Amcom)
Hillrom / Baxter	Stryker
	Vocera

Networking Infrastructure

Aruba	Extreme Networks
Cisco Spaces	Juniper Mist

Middleware / Integration Frameworks

IronBridge
Redox
Rhapsody

Workflow / Analytics

Dados Health	Olinqua
Helen	

Identity & Access Management

Active Directory
Single Sign-On (SSO)

HR

Dayforce	Workday
UKG	

Contact Us: www.centrak.com | marketing@centrak.com | 800-515-2928